



IRALINE Linear Hyaluronic Acid from **LKC Pharma** is made with only the best raw material.

**Hyaluronic Acid** is naturally synthesised by synovial cells.

In addition to being responsible for the viscoelastic properties of synovial fluid, it helps lubrication mechanisms under loading conditions, and has a protective function.

In degenerative joint disease, such a chronic primary osteoarthritis, synovial fluid has a reduced viscosity, decreasing its ability to lubricate the joint and absorb shocks.

Viscosupplementation with hyaluronic acid improves mobility and joint function, thereby reducing mechanical joint overload and, consequently, painful symptoms.



[www.lkcpharma.com](http://www.lkcpharma.com)

Via del Lavoro 4a/6  
20865 Usmate Velate (MB) - Italy  
+39 039.6317433

LKC Pharma injectable solutions are Class III Medical Devices (CE 0373). For medical use only.

  
Made in Italy



# IRALINE


## LINEAR HYALURONIC ACID



Iraline is a sterile, apyrogenic, resorbable solution made with hyaluronic acid created from non-animal sourced bacterial bio-fermentation.  
It is a gel made by linear hyaluronic acid.



## Suggested Treatments



## Tecnical Specifications\*

|                               |                                                                                     |
|-------------------------------|-------------------------------------------------------------------------------------|
| Composition                   | Sodium Hyaluronate, sterile saline solution                                         |
| HA Concentration              | 20mg/ml                                                                             |
| Raw Material molecular weight | 2,400 - 4,000 kDa                                                                   |
| Packaging                     | 1x2 ml                                                                              |
| Blister                       | 1 pre-filled, graduated, disposable.<br>Sterile glass syringe with Luer-Lok adapter |
| Suggested Application Depth   | Joint cavity                                                                        |
| Dosage                        | Once a week for 3-5 weeks                                                           |

\* Please refer to the package leaflet for full specifications